

TECHNOVENT

Technical/User Data Sheet **Procolour Sealant (P799)**

Overview

The Procolour extrinsic sealant is a dispersion of silicone rubber in a thin silicone based solvent. It is intended for sealing a silicone surface after painting (with Procolour and Extrinsic liquid paints).

The material is supplied in 2 x 20g bottles labeled part A and part B.

User Instructions

Dispense equal amounts of A and B into the supplied measuring cup and snap on lid. Shake vigorously to ensure thorough mixing and then paint onto the coloured silicone thinly. Alternatively weigh out equal parts of A and B and mix thoroughly.

After approx 20–30 minutes (depending on temperature of room) the sealant will become tacky – stipple the sealant surface with gauze in order to create a matt effect.

The sealant will be fully cured 30 minutes later – cure time can be speed up by heating in an oven.

TECHNOVENT

Technical/User Data Sheet Procolour Sealant (P799)

1. Identification of Substance & Company:

Product Name(s)	Silicone dispersion
Product Code(s)	P799
Medical Device Class	Class I (according to Annex IX of EC Directive 93/42/EEC)
Product Category/Classification	G5 (according to Annex IX of EC Directive 93/42/EEC)
Main Use	Sealing of painted silicones surfaces
Company	Technovent Ltd Unit 5 Kingsway Buildings, Bridgend South Wales, UK. CF31 3YH
Tel (& Emergency Tel)	(44) (0) 1656 768566
Fax	(44) (0) 1633 896588
Document Version	2
Print Date	12/09/2013



2. Composition / Information on ingredients:

Chemical Nature	Dispersion of platinum cure silicone rubber
Hazardous Components	Hexamethyl di siloxane

3. Hazards Identification:

Critical Hazards	FLAMMABLE
------------------	-----------



TECHNOVENT

Technical/User Data Sheet Procolour Sealant (P799)

4. First Aid Measures:

Eye Contact	May cause temporary irritation Irrigate with water for 15 mins, holding eyelid open
Skin Contact	Remove excess and wash with soap and water
Inhalation of Vapour	N/a
Ingestion	May result in gastric disturbances Seek medical advise

5. Fire Fighting Measures:

Hazards during firefighting	None
-----------------------------	------

6. Accidental Release Measures:

Personal Protection	Refer to Section 8
Environmental Precautions	None
Methods for clean-up	Collect with tissue or likewise and dispose as Section 13

7. Handling & Storage:

Handling	No smoking – do not breathe vapour Observe good housekeeping practises (HSE guidance note CS17)
Storage	Store at room temperature or refrigerated Keep away from sources of ignition in well ventilated plac

8. Exposure Controls / Personal Protection:

General Protective Measures	None deemed necessary
Components with workplace parameters	None, according to EU Directive 88/379/EC
Respiratory Protection	Not deemed necessary when used in small quantities
Hand Protection	Not deemed necessary, gloves recommended
Eye Protection	Not deemed necessary, eye glasses recommended
Skin Protection	Not deemed necessary
Hygiene Measures	Do not eat, drink or smoke when handling product

TECHNOVENT

Technical/User Data Sheet Procolour Sealant (P799)

15. Regulatory Information:

Symbol

FLAMMABLE



16. Other Information:

All information given in this Health & Safety Data Sheet is to the best of our knowledge true and accurate and is provided solely for making safety assessments. It is not a sales specification or an indication of suitability for a particular use nor does it replace the need for your own risk assessment.

All information presented in accordance with EC Directive 2001/58/EC and is correct at date of publication and is given in good faith but without warranty.

We cannot accept liability for any loss, damage or patent infringement resulting from the use of this information.

As with all materials, care should be exercised when handling.